

Medtronic

Engineering the extraordinary




05. Advanced automation features

MiniMed™ 780G system automatically adjusts* and corrects glucose levels every 5 minutes, as needed

 SmartGuard™ technology helps prevent highs and lows^{1,2,3}

Auto basal + Auto corrections



- A** Selection between a basal target of 5.5mmol/L (default setting), 6.1 mmol/L or 6.7mmol/L
- B** The auto correction target is set at 6.7 mmol/L
- C** Basal insulin adjusts every 5 minutes based on SG values
- D** Auto corrections delivered if the algorithm determines they are needed, and can be as frequent as every 5 minutes:
 - Max basal reached
 - $SG \geq 6.7\text{mmol/L}$No automatic corrections if Temp Target is set

* Refer to System User Guide - SmartGuard™ feature. Some user interaction required.

1. Carlson AL, et al. Diabetes Technol Ther. 2022;24(3):178-189.
2. Arrieta A, et al. Diabetes Obes Metab. 2022;24(7):1370-1379.
3. Collyns OJ, et al. Diabetes Care. 2021;44(4):969-975.

Background

Meal boluses are the main reason for post-prandial hypoglycaemia. The higher the insulin dose, the more risk of hypoglycaemia - even with a large meal¹

The MiniMed™ 780G system is outfitted with several modules for safety and efficacy.

- Safety of the MiniMed™ 780G algorithm for prevention of post meal and post correction bolus hypoglycaemia with [safe meal bolus](#) and the [safe correction bolus](#).
- Efficiency of the algorithm in reducing post-meal hyperglycaemia with the [meal detection](#) feature which is built into the safe correction bolus algorithm.




Safe meal bolus module


Safe meal bolus is a safety feature in the MiniMed™ 780G system that will reduce the amount of meal bolus insulin if the system predicts a hypoglycaemia would occur if the full bolus was given.

- Each time carbohydrates are entered into the SmartGuard™ bolus feature, the meal bolus is evaluated by the safe meal bolus feature.
- The safe meal bolus is a dynamic process to determine the meal bolus.

●..... The SmartGuard™ safe meal bolus feature is using:●




Safe bolus modules- brief overview



Safe meal bolus

Why is this occurring



- The safe meal bolus provides extra protection against hypoglycaemia
- Users may not be aware of auto-correction insulin or falling SG rate of change when delivering a meal bolus.
- Safe meal bolus is important because you may have more insulin in your system prior to meals as a result of automatic correction boluses.
- If the bolus isn't enough for the meal and glucose begins to rise, the pump can deliver automated insulin.
- The algorithm is programmed to prevent the need for preventive or “defensive” eating/carbohydrate ingestion.


Safe meal bolus

How often is this occurring

Rare occurrence?

Beside the common reasons for glucose rising post-meal (mistimed, miscalculated etc.), it can happen that the post-prandial glucose rise is caused by a bolus reduction of the safe meal bolus

- **Based on real-world data³, in less than 2% the bolus is reduced by more than 25%.**





< 2%
by more
than 25%

Safe meal bolus

Post prandial glucose rise?

1. Encourage patients to avoid focusing on single meals (most experience higher TIR, less TBR than previous therapy). Look at the “big picture” and the avoidance of the need for defensive eating.
2. Check the pump settings (AIT 2 hours, Target 5.5 mmol/L) and assure the max bolus is adequately set.
3. Assess bolus timing.
4. Optimise ICR as this is most impactful. (Ensure max bolus allows for ample bolus).
5. Ensure carbohydrate are entered and/or consume fewer carbohydrate. Generous carb estimates, individualised for patients, can assist once ICR is maximised according to the patient's experience.
6. Consider reduction in carb content and increase in proteins for recurrent issue.
7. If still no improvement, enter generous carbohydrate of 20-30% extra carbohydrate and adjust from there.

Case: Glucose rising post-meal





Most common reasons

- Insulin-to-carb ratio (ICR) not adequate?
- Bolus timing?
 - Check if there is already a pre-meal glucose rise
→ consider earlier timing of bolus
- Carb counting issues?
- Bolus omitted?

Recommended optimal settings* for better outcomes¹

In the real-world evidence

	All users	Recommended optimal setting users*
Number of users, n	61,481	4,392
Time in SmartGuard™, %	90.0	94.1
Mean SG, mmol/L	8.4	7.8
GMI**, %	6.9	6.7



1. Medtronic data on file: MiniMed™ 780G data uploaded voluntarily by 61,481 users in EMEA to CareLink™ Personal, from August 2020 to dec 2022.
*Users with recommended optimal settings are the ones using the combination of glucose target setting at 100 mg/dL and Active Insulin Time (AIT) at 2 hours for at least 95% of the time.

Safe meal bolus

Takeaway messages

- The algorithm is designed to keep patients safe by protecting against hypoglycaemia.
- The safe meal bolus will assist in reducing hypoglycaemia.
- Safe meal bolus is a dynamic process to determine the meal bolus.
- Trust the algorithm as it considers multiple factors. However, glucose spikes may occur but try not to focus on single meals.
- Focus on the big picture - more TIR, less TBR.
- Auto corrections help to lower glucose back into range.




Safe meal bolus

Takeaway messages

- To further mitigate post meal highs, assess the pump setting, try optimising ICR.
- Max bolus may need to be increased to allow for larger bolus.
- If the safe meal bolus feature is reducing the bolus size and post meal highs result, ensure carbohydrate are entered and/or consume fewer carbohydrate.
- Generous carb estimates, individualised for patients, can assist once ICR is maximised according to the patient's experience. For example, start by increasing carb count by 20-30%.




Example of Safe Meal Bolus



UK-IPT-2200173 | 2024

Medtronic

Example of Safe Meal Bolus



SmartGuard™ feature

Bolus timing – General recommendations


- Bolusing before meal is key to maximising time in SmartGuard™ Feature
- Deliver your bolus 10-20 minutes before eating, whether in SmartGuard™ Feature or Manual Mode

What can happen if a bolus is delivered during or after your meal?

- Glucose will begin to rise before the bolus is given
- Auto Basal will begin to increase, and Auto Corrections will be delivered, in an attempt to return glucose to your target
- Bolusing during or after meal can increase the risk for
 - high glucose values after the meal and/or
 - high glucose values followed by a precipitous drop in glucose

SmartGuard™ feature

Bolus timing: Real life example



- User eats before bolusing
 - Result: A sharp rise in sensor glucose (SG)
 - Auto Basal increase
- Carbs entered & bolus given mid to post - meal /snack
 - Result: A sharp drop in SG
 - SG often goes low to near-low "YO-YO" effect

SmartGuard™ feature

Bolus timing in other situations

Underestimating carbs

- Enter and bolus for additional carbs at the end of the meal.
- If > 1 hour since eating, enter current BG and give the system recommended correction bolus

Unsure of how much will be eaten

- Enter and bolus pre-meal for the grams that you certain will eat
- If additional grams are consumed, enter carbs and bolus for additional grams as eaten

Missing a bolus

- If ≤ 1 hour after eating, enter half the grams eaten and recheck BG in an hour
 - If glucose is still high, enter current BG and give the system recommended correction bolus
- If > 1 hour after eating, enter the BG and give the system recommended correction bolus.

Meals high in fat and protein

How to bolus when in SmartGuard™ feature

- With usual pump therapy or in Manual Mode with the MiniMed™ 780G System you are able to use bolus options to match the insulin timing with the effect of the meal on your glucose level
- When in SmartGuard™ Feature, the algorithm adjusts in real-time to compensate for small miscalculations in carbs, high-fat meals, and delayed digestion.
- You may need less insulin for meals rich in fat and protein because of the slower digestion of carbs
- The MiniMed™ 780G System will use the Insulin-Carb-Ratio (ICR) from your Bolus Wizard™ settings to calculate the bolus
 - You only need to enter BG and carbs
- Keep in mind, dosing adjustments are conservative and not intended to correct for large mismatches in carbs or a missed bolus



Basal delivery is adjusted based on rise or fall of glucose.

SmartGuard™ feature

Bolusing for high fat / high protein meals


Please note:

- Auto Basal can usually accommodate glucose fluctuations associated with high fat / high protein meals
- Dual Wave™ and Square Wave™ boluses are not available in SmartGuard™ Feature

If high glucose values do occur:

- Suggest entering 50%-75% of grams before eating and entering the remainder 1-2 hours after the meal
- Depending on meal composition, some people may require the majority of carbs to be entered before eating, with additional grams added later to compensate for the high fat content and slow digestion of the meal

MiniMed™ 780G system with SmartGuard™ feature helps to simplify your patients' exercise¹



1. Lee et al. Diabetes Care. Volume 43, February 2020 480-483
This representation provides guidelines. Consult with your healthcare professional on how to manage exercise with your device.

Low glucose

How to treat when using SmartGuard™



Smaller amounts of carbohydrate are generally needed due to the reduction or suspension of insulin delivery before the onset of hypoglycaemia.³
Half of your usual treatment is usually a good starting point.⁴

Check
glucose
in 15
minutes

Repeat
treatment
if glucose still
below
3.9 mmol/L



Information contained herein does not replace the recommendations of your healthcare professional.



Simplified Meal Management


Delivers proven clinical outcomes in a wide T1D population*



Improvements in TIR were significantly **higher with flex** compared to simplified meal announcement, but **HbA1c and mean SG** were **comparable**



Fix group reached international **targets of glycaemic control**



The AHCL algorithm is capable of overcoming **lower precision** in carb counting



Meal management with precise carbohydrate counting further improves outcomes, **carbohydrate estimation skills remain important**

*The MiniMed™ 780G insulin pump is indicated for use by patients aged 7-80 years with Type 1 diabetes

1. Petrovski G, et al. Diabetes Care 2023;46(3):544-550

Q&A

References

1. Slattery D, et al. Optimal prandial timing of bolus insulin in diabetes management: a review. Diabet. Med. 2018; 35: 306-316
2. Castañeda J, Mathieu C, Aanstoot HJ, Arrieta A, Da Silva J, Shin J, Cohen O. Predictors of time in target glucose range in real-world users of the MiniMed 780G system. Diabetes Obes Metab. 2022 Nov;24(11):2212-2221. doi: 10.1111/dom.14807. Epub 2022 Aug 1. PMID: 35791621
3. MiniMed™ 780G data uploaded voluntarily by 2944 users in EMEA to CareLink™ Personal, from 27 August 2020 to 11 January 2021

Disclaimer

MiniMed™ 780G system

This material does not replace or supersede the instructions for use. It should not be considered the exclusive source of information and should be used in conjunction with the User Guide.

See the User Guide for detailed information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic website at [Medtronic diabetes](https://www.medtronic.com/diabetes)

UK-IPT-2300083 © 2023 Medtronic. All rights reserved. Medtronic and the Medtronic logo are trademarks of Medtronic.

TM* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

ACCU-CHEK, ACCU-CHEK GUIDE, ACCU-CHEK FASTCLIX, and ACCU-CHEK SOFTCLIX are trademarks of Roche Diabetes Care.

DreaMed Diabetes is a trademark of DreaMed Diabetes, Ltd. The MiniMedTM 780G system algorithm includes technology developed by DreaMed Diabetes